

## NeuroTrac Labour TENS

### Safe and Drug Free Pain Relief in Childbirth



The **Neurotrac Labour Tens** is an advanced digital dual channel TENS machine that features a dedicated Labour Pain programme. Three different treatment modes provide maximum coverage of the labour pains spectrum. The NeuroTrac Labour Tens can help to relieve obstetric pain and all kinds of aches associated with the childbirth. 80% efficacy rate has been observed.



**Neurotrac Labour Tens presents an effective alternative to pharmacological (drug based analgesia) pain relieving therapies. NeuroTrac Labour Tens is a viable alternative for all mums to be, who wish to deliver their babies naturally and as comfortably as possible and without unwanted side effects.**

The unit is controlled by a hand held push button control that allows easy and convenient change of stimulation modes for the best pain control during labour.

NeuroTrac Labour Tens uses either four long rectangular self adhesive skin electrodes that are 100 x 50 cm in size or two very long self-adhesive rectangular electrodes, which measure 204 x 50 cm.



**The four electrodes (100 x 50 cm) should be placed between the 10th thoracic vertebrae and the 1st lumbar vertebrae as pictures on the left.**

Alternatively:



**The two long electrodes (204 x 50 cm) should be placed between the 10th thoracic vertebrae and the 2nd sacral vertebrae as pictures on the right.**

Clinical evidence suggests that the optimum pain control is achieved, when the stimulator is used from the first onset of the labour pain and for two to three hours after the birth. The low frequency stimulation is to be used during the periods between contractions, whilst the higher frequency should be reserved for the times of the strong contractions.

#### NeuroTrac Labour Tens offers the following specifications:

- Easy to use digital dual channel unit for obstetric pain relief during labour
- 1 dedicated labour pain programme for ease of use
- Remote boost switch for additional pain control at times of intense pain
- Continuous mode (Boost): 90 Hz at 220  $\mu$ S
- Burst mode (Rest): 2 Hz Bursts. Bursts of 9 pulses at 200  $\mu$ S and 150 Hz repeating twice per second
- Constant current with extended battery life
- Automatic switch off after 4 minutes, if the unit is not in use



Win Health Ltd., Brockhirst, Oxnam Road, Jedburgh, TD8 6QN  
Tel. 01835 864866; Fax: 01835 863238  
sales@winhealth.co.uk  
[www.winhealth.co.uk](http://www.winhealth.co.uk)

## Neurotrac Labour TENS



*Please note that the NeuroTrac Labour Tens has been designed to help women during labour and therefore has only a limited use for other pain relief purposes.*



**NeuroTrac Labour TENS is a CE marked Medical Device.**

**Neurotrac Labour TENS is presented in a convenient carry pouch with 1 set of 4 skin electrodes (100 x 50 cm), 1 set of 2 skin electrodes (130 x 50 cm), a push button control unit, 2 sets of leads and a 9 v alkaline battery. Easy to follow user instructions manual included.**

**2 years warranty on materials, components and workmanship.**

**Womens Waterworks Curing Incontinence by Dr Pauline Chairelli that provides all Mums with an important information on pelvic floor health after the childbirth is included with each Neurotrac Labour TENS set.**

**NeuroTrac Labour Tens is affordably priced at £39.99 including VAT.**

*Please note that there is no VAT exemption available on this product! Please note that the NeuroTrac Labour Tens that we sell are brand new and NOT ex-rental! Additional/Spare skin electrodes are inexpensive and readily available.*



**Win Health Ltd.  
Brockhirst  
Oxnam Road  
Jedburgh, TD8 6QN  
Roxburghshire  
Tel. 01835 864866  
Fax: 01835 863238**

**[sales@winhealth.co.uk](mailto:sales@winhealth.co.uk)  
[www.winhealth.co.uk](http://www.winhealth.co.uk)**



Win Health Ltd., Brockhirst, Oxnam Road, Jedburgh, TD8 6QN  
Tel. 01835 864866; Fax: 01835 863238  
[sales@winhealth.co.uk](mailto:sales@winhealth.co.uk)  
[www.winhealth.co.uk](http://www.winhealth.co.uk)